

### VideometerLab for combined spectral reflectance and fluorescence imaging







## **ABOUT US**



- Spectral imaging company
- Founded 1999
- Products
  - Lab instruments,
  - Turn-key in-line systems, and
  - R&D projects
- App. 700 imaging R&D projects since 2000
- In-line 24/7 spectral imaging since 2002
- Based in Copenhagen, Denmark
- Partnerships worldwide

## **OUR LEGACY**





#### The beginnings

Videometer was co-founded by Jens Michael Carstensen and 7-Technologies in 1999, as a spinoff from the Technical University of Denmark. The first patent application was filed.

#### **Project-based**

In 2000, Videometer began its project-based activity. During these years, the company's main focus was set on custom-made vision systems for in-line and on-line quality control.

#### A new era

In 2018, Videometer's structure underwent new developments both in terms of strategy and structure. This year marked the beginning of a new era for the company, in terms of focus on instruments.

#### Today

Today, Videometer is a leading provider of spectral imaging solutions worldwide, selling both spectral imaging instruments and custom-made vision systems. Videometer is synonym of excellence and innovation in its field.





## **SPECTRAL IMAGING**



#### WHAT COLOR IS THE CAR?





APPEARANCE CHEMISTRY Х PHYSICS Х ENVIRONMENT Х ILLUMINATION

#### LED BAND SEQUENTIAL SPECTRAL IMAGING





Camera and lens

Emission filter changer

Integrating sphere

LEDs of multiple wavelengths

Sample is placed in target opening

Backlight or background



- LEDs: Stable, durable, large selection, rapidly developing technology
- Up to 20 different high-resolution bands acquired sequentially in 0.5-1.0 seconds
- May be combined with emission filters, backlight, and darkfield illuminant
- Combined reflectance spectral imaging and fluorescence spectral imaging possible!

## SPECTRAL IMAGE







N images obtained at N wavelengths

Microbial and plant as metabolites	Accurate color assessment and pigment concentration	Pigment baseline, moisture, fat, etc.		Spectral image is typically a large data structure of 100 MB to 10 GB
--	---	--	--	---

### MULTISPECTRAL IMAGING





## **EMISSION FILTER CHANGER**



- 5 emission filter apertures where the first is usually left empty for reflectance imaging
- Typically 4 apertures for long pass (LP) filters (fluorescence) or short pass (SP) filters (upconversion). Band pass filters (BP may also be used)
- All filters must be 25 mm diameter and up to 3 mm thickness. Possible suppliers are e.g.
  - Andover Corporation (thickness is typically 1.5 mm unmounted)
    <a href="http://www.andovercorp.com/Web">http://www.andovercorp.com/Web</a> store/Edge Filters/steep edge filters.php
  - Semrock (be aware that mounted filters are 3.5 mm thick) <a href="http://www.semrock.com/filtersRefined.aspx?id=21&page=1&so=0&recs=100">http://www.semrock.com/filtersRefined.aspx?id=21&page=1&so=0&recs=100</a>
  - Edmund Optics (thickness is typically 3 mm unmounted)

http://www.edmundoptics.com/optics/optical-filters/longpass-edge-filters/highperformance-od-4-longpass-filters/3044

## **EMISSION FILTER SET**



- Standard emission filter sets
  - LP OD4 400 nm, 500 nm, 600 nm, 700 nm
  - LP OD4 450 nm, 550 nm, 650 nm, 700 nm
- Customized filter sets available
- Easy filter set configuration in software
- Adaptive setting of exposure time and optimal dynamic range
- One image obtained from each combination of LED excitation strobe and emission filter. LED strobes in the pass band of emission filters are usually omitted.



#### WHEAT WITH BUNT 4 SPORED PER SEED



### **FUSARIUM FLUORESCENCE**





We measure what you see – and beyond



### **ACROSPIRE IN BARLEY MALT**



Acrospire visible

We measure what you see – and beyond

## **OUR VALUES**





Responsible Consumption and Production



Good Health and Well-Being



Life Below Water



Decent Work and Economic Growth



Partnership for the Goals



# **THANK YOU!**

